

# K300 PS Jr



**\* Indicates a Required Field**

## Dealer Information



\* **Contact:** \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* **Dealer Name:** \_\_\_\_\_

Address: \_\_\_\_\_

\* **City:** \_\_\_\_\_

\* **State/Zip:** \_\_\_\_\_

\* **Phone#:** \_\_\_\_\_

\* **Fax#:** \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* **First Name:** \_\_\_\_\_

\* **Last Name:** \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* **Weight:** \_\_\_\_\_

\* **Height:** \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\* **I) Knee to Foot:** \_\_\_\_\_



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobilus.com](mailto:sales@permobilus.com)

**Permobil Inc.**  
6961 Eastgate Blvd.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.com](http://www.permobil.com)

Prices effective **August 9, 2010.**

# K300 PS Jr

## Chair Model and Colors

### CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103573-99-0	K300 Base PS - C300 PS0 <i>Price Includes: FWD C300 Power Base Set Up For No Power Seat Functions, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,845.00	K0848
<input type="checkbox"/> I10431	K300 Base PS - C300 PS1 <i>Price Includes: FWD C300 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,845.00	K0856
<input type="checkbox"/> I10520	K300 Base PS - C300 PS2 <i>Price Includes: FWD C300 Power Base Set Up For Two Power Seat Functions, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,845.00	K0861

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge	
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103576-99-0	C300/K300 Shroud - Onyx Black	No Charge	
<input type="checkbox"/> I101383-99-0	C300/K300 Shroud - Crystal Blue	No Charge	
<input type="checkbox"/> I103577-99-0	C300/K300 Shroud - Platinum Silver	No Charge	
<input type="checkbox"/> I101385-99-0	C300/K300 Shroud - Green Metallic	No Charge	
<input type="checkbox"/> I103575-99-0	C300/K300 Shroud - Titan Blue	No Charge	
<input type="checkbox"/> I103578-99-0	C300/K300 Shroud - Pop Star Pink	No Charge	
<input type="checkbox"/> I103574-99-0	C300/K300 Shroud - Chili Pepper Red	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge	
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge	

## Base Options

### MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101995-99-0	Fixed Seat Tube for Tilt	No Charge	
<input type="checkbox"/> I103779-99-0	Power Adjustable Seat Height - 8" Travel	2,445.00	E2300

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IM34 SLDG	Batteries - Group 34, 12V 60 Amp Hours, Installed	785.00	K0108
<input type="checkbox"/> I102768-99-0	Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent &amp; a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics).</i>	402.00	E1029
<input type="checkbox"/> I0036-7002A	Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	1,170.00	K0108
<input type="checkbox"/> I103564-99-0	Light Kit for C300 w/ R-net Electronics <i>Shrouds with lights cannot be ordered with VR2 Electronics.</i>	972.00	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	443.00	

# K300 PS Jr

## Permolock

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to:</i> NAME: _____ ADDR: _____	1,806.00	
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	605.00	

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> I103544-99-0	R-net Remote Joystick w/ Mono Jacks	539.00	
<input type="checkbox"/> ID51040	R-net Remote Joystick w/ Toggle Switches	616.00	
<input type="checkbox"/> ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle	743.00	
<input type="checkbox"/> ID51153	R-net Remote Color Screen Joystick w/ Mono Jacks	850.00	
<input type="checkbox"/> ID51149	R-net Remote Color Screen JS w/ Toggle Switches	890.00	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Speed Pot & Single Toggle	890.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode &amp; on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	616.00	E2373
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	198.00	
<input type="checkbox"/> ID51108	R-net LED Remote Joystick Module <i>Standard Joystick Module. This Joystick does NOT have an LCD Display and does NOT contain mono jacks. Limited to five profiles with this Joystick Module. Please note if ordering this item, the ICS Switchbox must be mounted onto either the armrest pad or the backrest. Cannot order PN: I1822845 or I103469-99-0 with this style of joystick.</i>	No Charge	

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	96.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	141.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	96.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	129.00	E2323
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	129.00	E2323

### CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	297.00	E1028

## K300 PS Jr

### Joystick Options

#### CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	297.00	E1028
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	495.00	E1028
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	495.00	E1028

#### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0) OR when ordering Alternative Drive Packages.</i>	545.00	E2377
<input checked="" type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	475.00	E2313
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i>	1,103.00	E2331
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	809.00	
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	1,012.00	E2351

### Seating System

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103712-99-0	PS Junior Seat Frame - Only <i>Weight limit 165lbs.</i>	No Charge	
<input type="checkbox"/> I103715-99-0	PS Jr Seat w/ 45° Manual Tilt <i>Must be ordered with Stroller Handles (I103811-99-0). Weight limit 165lbs.</i>	625.00	K0108
<input type="checkbox"/> I103716-99-0	PS Junior Seat Frame with 45° Power Tilt <i>Weight limit 165lbs.</i>	5,495.00	E1002

### Seating System Options

#### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	1,495.00	E2310
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	2,925.00	E2311

#### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	249.00	

# K300 PS Jr

## Seating System Options

### CHOOSE MOUNTING POSITION OF SEAT FUNCTION CONTROL KIT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	127.00	

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,002.00	K0108
<input type="checkbox"/> I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	1,002.00	

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Solid Back, Sling Back, or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10448	PS Jr Curved Metal Backrest <i>Backrest width will match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the backrest.</i>	No Charge	
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge	

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10432	PS Jr Seat Width 11"	No Charge	
<input type="checkbox"/> I10433	PS Jr Seat Width 12"	No Charge	
<input type="checkbox"/> I10434	PS Jr Seat Width 13"	No Charge	
<input type="checkbox"/> I10435	PS Jr Seat Width 14"	No Charge	
<input type="checkbox"/> I10436	PS Jr Seat Width 15"	No Charge	
<input type="checkbox"/> I10437	PS Jr Seat Width 16"	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10438	PS Jr Seat Depth 10"	No Charge	
<input type="checkbox"/> I10439	PS Jr Seat Depth 11"	No Charge	
<input type="checkbox"/> I10440	PS Jr Seat Depth 12"	No Charge	
<input type="checkbox"/> I10441	PS Jr Seat Depth 13"	No Charge	
<input type="checkbox"/> I10442	PS Jr Seat Depth 14"	No Charge	
<input type="checkbox"/> I10443	PS Jr Seat Depth 15"	No Charge	
<input type="checkbox"/> I10444	PS Jr Seat Depth 16"	No Charge	
<input type="checkbox"/> I10445	PS Jr Seat Depth 17"	No Charge	

## K300 PS Jr

### Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10446	PS Jr Seat Depth 18"	No Charge	

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10038	Ergonomic Seat Cushion Leatherette <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	483.00	E2601,E2602
<input type="checkbox"/> I10276	Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Not available with seat widths or depths greater than 14".</i>	483.00	

### Armrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104648-99-0	Height Adj Arm Assembly - PS Jr w/ UT Arm Pads <i>Newly Redesigned! Includes: 40° of angle adjustability, unlimited rotational adjustment, and 2" of offset adjustment medially and laterally. Height adjustment in 1" increments from 7" - 11" from the seat pan. Includes "C" Sized UT Armrest Pad (8"L x 3.5"W). Height Adjustable Armrest Assembly is no charge if ordering Power Tilt. Other sized armpads are be available upon request.</i>	289.00	E0973
<input type="checkbox"/> I10519	Height Adj Arm Assembly - PS Jr w/ UT Arm Pads - NC <i>No Charge when ordering Powered Tilt.</i>	No Charge	

### Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10586	70° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate.</i>	No Charge	
<input type="checkbox"/> I10587	90° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate.</i>	No Charge	
<input type="checkbox"/> I10625	Aeron Legrest Assembly for PS Jr Seat <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes Narrow, Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	300.00	

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104684-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 &amp; I10587).</i>	No Charge	
<input type="checkbox"/> I104686-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 &amp; I10587).</i>	No Charge	
<input type="checkbox"/> IFT202SR	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	229.00	K0108

### Accessories

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge	
<input type="checkbox"/> I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	115.00	K0108
<input type="checkbox"/> I10623	Bodypoint Padded Hip Belt UniTrack - Extra Small	115.00	
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	239.00	K0108
<input type="checkbox"/> I100749-99-0	Miniflex Headrest Removable - Leatherette <i>Must also select PN I10552.</i>	268.00	E0955

# K300 PS Jr

## Accessories

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	205.00	K0108
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	431.00	E1028
<input type="checkbox"/> I103811-99-0	Stroller Handle for PS Jr Seating System <i>Angle adjustable. This item Must be selected when ordering manual tilt (I103715-99-0).</i>	185.00	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	229.00	
<input type="checkbox"/> I101473-99-0	Medical Necessity Bag Clips - Carabiner Style	76.00	
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	199.00	E0960
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	242.00	E0950
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	301.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I103798-99-0	After Market Seating System Adapter <i>Used to connect after market seating systems to the UniTrack Accessory Rail on the PS Jr Seat Frame. Includes Mounting Hardware.</i>	76.00	K0108

Custom Swing Away Laterals Available Upon Request

## Alternative Drive Systems

### (1) CHOOSE ALTERNATIVE DRIVE INPUT:

Part Number	Description	Price	HCPCS Code
	R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Adjustable/Removable Alt Drive Hardware, AND (3) Omni or Input Module.</i>		
<input type="checkbox"/> I10558	Gimble Joystick set up for Armrest Mount	1,412.00	E2373
<input type="checkbox"/> I10559	Gimble Joystick set up for Midline Mount	1,412.00	E2373
<input type="checkbox"/> I10560	Gimble Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	1,412.00	E2373
<input type="checkbox"/> I10561	Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	1,582.00	E2373
<input type="checkbox"/> I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,582.00	E2373
<input type="checkbox"/> I10564	Mini Joystick set up for Armrest Mount	2,685.00	E2312
<input type="checkbox"/> I10565	Mini Joystick set up for Midline Mount	2,685.00	E2312
<input type="checkbox"/> I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,745.00	E2312
<input type="checkbox"/> I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,745.00	E2312
<input type="checkbox"/> I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,265.00	E2312

# K300 PS Jr

## Alternative Drive Systems

### (1) CHOOSE ALTERNATIVE DRIVE INPUT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,375.00	E2328
<input type="checkbox"/> I10570	Head Array - Adult 5 Switch HDR-A5	4,375.00	E2330
<input type="checkbox"/> I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	4,175.00	E2330
<input type="checkbox"/> I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	4,375.00	E2330
<input type="checkbox"/> I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,475.00	E2325
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,963.00	E2325,E2326

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10076	Right Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10077	Left Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10074	Right Side Mount For Omni	No Charge	
<input type="checkbox"/> I10075	Left Side Mount For Omni	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Panel on Right	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Panel on Left	No Charge	
<input type="checkbox"/> I10080	Mount Seat Control Panel on Backrest	No Charge	

## Alt Drive Systems Setup

### (2) CHOOSE ADJ ALT DRIVE HARDWARE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10577	Adjustable Removable Hardware - Alt Drive	216.00	E1028

### (3) CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	1,496.00	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,012.00	E2351

## Alt Drive Accessories

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	151.00	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	113.00	K0108
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	110.00	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	524.00	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	603.00	
<input type="checkbox"/> I5018	Wafer Board w/5 Switches	506.00	

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*