

# K300 PS Jr



## Dealer Information

Contact: \_\_\_\_\_  
Dealer Code: \_\_\_\_\_  
Dealer Name: \_\_\_\_\_  
Address: \_\_\_\_\_  
City: \_\_\_\_\_  
Province/Postal Code: \_\_\_\_\_  
Phone#: \_\_\_\_\_  
Fax#: \_\_\_\_\_  
PO#: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*


First Name: \_\_\_\_\_  
Last Name: \_\_\_\_\_  
Diagnosis: \_\_\_\_\_  
Age: \_\_\_\_\_

## Client Measurements

**Weight:** \_\_\_\_\_

**Height:** \_\_\_\_\_

- A) Top of Shoulders: \_\_\_\_\_  
B) Chest Depth: \_\_\_\_\_  
C) Chest Width: \_\_\_\_\_  
D) Seat Depth: \_\_\_\_\_  
E) Top of Head: \_\_\_\_\_  
F) Elbow to Hand: \_\_\_\_\_  
G) Seat Pan to Elbow: \_\_\_\_\_  
H) Hip Width: \_\_\_\_\_  
I) Knee to Foot: \_\_\_\_\_



Display ADP Codes on order/quote: \_\_\_\_\_

Please send order/quote to fax#: **(800) 231-3256**

Canadian prices effective **October 28, 2008** and payable in Canadian dollars. All specifications and prices are subject to change without notice. **Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.**

**Permobil Ltd.**  
Toronto Canada  
C/O Permobil Inc  
6961 Eastgate Blvd.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
www.permobil.com

**Base***Select Options:*

Part Number	Description	MSRP - CAD	ADP Code
I10431	K300 Base PS - C300 PS1 <i>Price Includes: FWD C300 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,995.00 <input type="checkbox"/>	WPK31P010

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge <input type="checkbox"/>	
I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge <input type="checkbox"/>	

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I103576-99-0	C300/K300 Shroud - Onyx Black	No Charge <input type="checkbox"/>	
I101383-99-0	C300/K300 Shroud - Crystal Blue	No Charge <input type="checkbox"/>	
I103577-99-0	C300/K300 Shroud - Platinum Silver	No Charge <input type="checkbox"/>	
I101385-99-0	C300/K300 Shroud - Green Metallic	No Charge <input type="checkbox"/>	
I103575-99-0	C300/K300 Shroud - Titan Blue	No Charge <input type="checkbox"/>	
I103578-99-0	C300/K300 Shroud - Pop Star Pink	No Charge <input type="checkbox"/>	
I103574-99-0	C300/K300 Shroud - Chili Pepper Red	No Charge <input type="checkbox"/>	

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge <input type="checkbox"/>	
I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge <input type="checkbox"/>	

**Base Options***Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I101995-99-0	Fixed Seat Tube for Tilt	No Charge <input type="checkbox"/>	
I103779-99-0	Power Adjustable Seat Height - 8" Travel <i>To assist in transfers and improve access for function.</i>	2,495.00 <input type="checkbox"/>	

*Select Options:*

Part Number	Description	MSRP - CAD	ADP Code
IM34 SLDG	Batteries - Group 34, 12V 60 Amp Hours, Installed	762.00 <input type="checkbox"/>	
I102768-99-0	Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent &amp; a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics).</i>	395.00 <input type="checkbox"/>	WEPY
I0036-7002A	Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0).</i>	1,250.00 <input type="checkbox"/>	WEPQ
I103564-99-0	Light Kit for C300 w/ R-net Electronics <i>Shrouds with lights cannot be ordered with VR2 Electronics.</i>	955.00 <input type="checkbox"/>	
ID51113	Programming Key for R-net - Permobil Version	435.00 <input type="checkbox"/>	

**Joystick Options***Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10045	Cone Shaped Joystick Knob - Standard.	No Charge <input type="checkbox"/>	

## Joystick Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I100960-99-0	Chin Cup For Joystick Knob	94.00 <input type="checkbox"/>	
I100961-99-0	"T" Handle For Joystick Knob	139.00 <input type="checkbox"/>	
I100962-99-0	"Large Ball" For Joystick Knob	94.00 <input type="checkbox"/>	
I100963-99-0	"Softball" For Joystick Knob	127.00 <input type="checkbox"/>	
I100964-99-0	"Mushroom" For Joystick Knob	127.00 <input type="checkbox"/>	
I100965-99-0	"Stick" For Joystick Knob	127.00 <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	325.00 <input type="checkbox"/>	WEPM
I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	325.00 <input type="checkbox"/>	WEPM
I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	375.00 <input type="checkbox"/>	WEPM
I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	375.00 <input type="checkbox"/>	WEPM
ID51108	R-net LED Remote Joystick Module <i>Standard Joystick Module. This Joystick does NOT have an LCD Display and does NOT contain mono jacks. Limited to five profiles with this Joystick Module. Must order expandable controller (PN I10411) with this part. Please note if ordering this item, the ICS Switchbox must be mounted onto either the armrest pad or the backrest. Cannot order PN: I1822845 or I103469-99-0 with this style of joystick.</i>	No Charge <input type="checkbox"/>	
I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. This is an expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode &amp; on/off.</i>	1,650.00 <input type="checkbox"/>	
I103544-99-0	R-net Remote Joystick w/ Mono Jacks <i>Expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	1,495.00 <input type="checkbox"/>	
ID51040	R-net Remote Joystick w/ Toggle Switches <i>Expandable/remote joystick with toggle switches. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	1,595.00 <input type="checkbox"/>	WEX1
ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle <i>Expandable/remote joystick with a single toggle switch for modes/profiles &amp; a speed pot. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	1,595.00 <input type="checkbox"/>	WEX1
I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i>	No Charge <input type="checkbox"/>	
ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 controllers.</i>	1,084.00 <input type="checkbox"/>	
I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	1,450.00 <input type="checkbox"/>	

## Seating System

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103712-99-0	PS Junior Seat Frame - Only <i>Weight limit 165lbs.</i>	No Charge <input type="checkbox"/>	

## Seating System

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103715-99-0	PS Jr Seat w/ 45° Manual Tilt	525.00 <input type="checkbox"/>	
I103716-99-0	PS Junior Seat Frame with 45° Power Tilt <i>Must be ordered with Stroller Handles (I103811-99-0). Weight limit 165lbs. Weight limit 165lbs.</i>	5,595.00 <input type="checkbox"/>	SEC1P0205

## Seating System Options

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10418	Single Seat Function Control Switch for R-net	No Charge <input type="checkbox"/>	
I10419	Single Seat Function Control Kit for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	1,378.00 <input type="checkbox"/>	
I10585	Multiple Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	2,295.00 <input type="checkbox"/>	SEC1P1005
I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge <input type="checkbox"/>	
I102799-99-0	Toggle Switches on Seat Function Control Kit	245.00 <input type="checkbox"/>	
I103469-99-0	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge <input type="checkbox"/>	
I103470-99-0	Mount ICS Switchbox on Backrest	125.00 <input type="checkbox"/>	
I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,084.00 <input type="checkbox"/>	

## Backrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable	No Charge <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10448	PS Jr Curved Metal Backrest <i>Backrest width will match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the backrest.</i>	No Charge <input type="checkbox"/>	
I10449	PS Jr - Omit Backrest	No Charge <input type="checkbox"/>	

## Seat Sizing

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10432	PS Jr Seat Width 11"	No Charge <input type="checkbox"/>	
I10433	PS Jr Seat Width 12"	No Charge <input type="checkbox"/>	
I10434	PS Jr Seat Width 13"	No Charge <input type="checkbox"/>	
I10435	PS Jr Seat Width 14"	No Charge <input type="checkbox"/>	
I10436	PS Jr Seat Width 15"	No Charge <input type="checkbox"/>	

**Seat Sizing***Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10437	PS Jr Seat Width 16"	No Charge <input type="checkbox"/>	

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10438	PS Jr Seat Depth 10"	No Charge <input type="checkbox"/>	
I10439	PS Jr Seat Depth 11"	No Charge <input type="checkbox"/>	
I10440	PS Jr Seat Depth 12"	No Charge <input type="checkbox"/>	
I10441	PS Jr Seat Depth 13"	No Charge <input type="checkbox"/>	
I10442	PS Jr Seat Depth 14"	No Charge <input type="checkbox"/>	
I10443	PS Jr Seat Depth 15"	No Charge <input type="checkbox"/>	
I10444	PS Jr Seat Depth 16"	No Charge <input type="checkbox"/>	
I10445	PS Jr Seat Depth 17"	No Charge <input type="checkbox"/>	
I10446	PS Jr Seat Depth 18"	No Charge <input type="checkbox"/>	

*Select Options:*

Part Number	Description	MSRP - CAD	ADP Code
I10038	Ergonomic Seat Cushion Leatherette <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	475.00 <input type="checkbox"/>	SES1P0005
I10276	Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Not available with seat widths or depths greater than 14".</i>	475.00 <input type="checkbox"/>	

**Armrest Options***Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10588	Height Adjustable Armrest Assembly for PS Jr - NC <i>Includes: 40° of angle adjustability, unlimited rotational adjustment, and 2" of offset adjustment medially and laterally. Height adjustment in 1" increments from 7" - 11" from the seat pan. Height Adjustable Armrest Assembly is no charge if ordering Power Tilt.</i>	No Charge <input type="checkbox"/>	
I10589	Height Adjustable Armrest Assembly for PS Jr <i>Charged item when ordering non-powered tilt systems.</i>	284.00 <input type="checkbox"/>	

**Legrest Options***Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10586	70 Degree Swing Away Legrests - PS Jr	No Charge <input type="checkbox"/>	
I10587	90 Degree Swing Away Legrests - PS Jr	No Charge <input type="checkbox"/>	

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10590	Multi-Axis Angle Adjust Footplates - 5"W x 8"D	No Charge <input type="checkbox"/>	

*Select Options:*

Part Number	Description	MSRP - CAD	ADP Code
-------------	-------------	------------	----------

**Accessories***Select Options:*

Part Number	Description	MSRP - CAD	ADP Code
I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge <input type="checkbox"/>	
I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	119.00 <input type="checkbox"/>	
I100519-99-0	Positioning Belt Retractable	235.00 <input type="checkbox"/>	

## Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I100749-99-0	Miniflex Headrest Removable - Leatherette <i>Must also select PN I10552.</i>	410.00 <input type="checkbox"/>	
IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	197.00 <input type="checkbox"/>	
I10552	Adj Removable - Headrest Hardware (1 Each)	110.00 <input type="checkbox"/>	SEHND0020
I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	410.00 <input type="checkbox"/>	SELND0015
I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	410.00 <input type="checkbox"/>	SELND0015
I10550	Adj Removable - Thigh Support Hardware (Pair)	110.00 <input type="checkbox"/>	
I103811-99-0	Stroller Handle for PS Jr Seating System <i>Angle adjustable. This item Must be selected when ordering manual tilt (I103715-99-0).</i>	182.00 <input type="checkbox"/>	
I101473-99-0	Back Pack Clips Carabiner Style	75.00 <input type="checkbox"/>	
I102059-99-0	Medical Necessities Bag <i>Includes Back Pack Clips (Carabiner Style).</i>	265.00 <input type="checkbox"/>	
I102697001331	Medical Necessities Bag - Case Logic	265.00 <input type="checkbox"/>	
I103798-99-0	After Market Seating System Adapter <i>Used to connect after market seating systems to the UniTrack Accessory Rail on the PS Jr Seat Frame. Includes Mounting Hardware.</i>	75.00 <input type="checkbox"/>	

## Alternative Drive Systems

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Expandable Controller, (3) Adjustable/Removable Alt Drive Hardware, AND (4) Omni or Input Module.</i>		
	(1) INPUT DEVICES FOR ALTERNATIVE DRIVE CONTROL		
I10558	Gimble Joystick set up for Armrest Mount	1,388.00 <input type="checkbox"/>	WEX2
I10559	Gimble Joystick set up for Midline Mount	1,388.00 <input type="checkbox"/>	WEX2
I10560	Gimble Joystick set up for Swing Away Chin Boom	1,388.00 <input type="checkbox"/>	WEX2
I10561	Gimble Joystick w/ Haims Harness - Small <i>Custom charge may be applied if ordering lower backrest heights.</i>	1,555.00 <input type="checkbox"/>	WEX2
	<i>For neck width up to 5".</i>		
I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,555.00 <input type="checkbox"/>	WEX2
I10563	Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00 <input type="checkbox"/>	WEX2
I10564	Mini Joystick set up for Armrest Mount	2,035.00 <input type="checkbox"/>	WEX2
I10565	Mini Joystick set up for Midline Mount	2,035.00 <input type="checkbox"/>	WEX2
I10566	Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	2,035.00 <input type="checkbox"/>	WEX2
I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,135.00 <input type="checkbox"/>	WEX2
I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,135.00 <input type="checkbox"/>	WEX2
I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00 <input type="checkbox"/>	

## Alternative Drive Systems

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,085.00 <input type="checkbox"/>	WEX3
I10570	Head Array - Adult 5 Switch HDR-A5	4,075.00 <input type="checkbox"/>	WEX4
I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	4,075.00 <input type="checkbox"/>	WEX4
I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	4,075.00 <input type="checkbox"/>	WEX4
I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,409.00 <input type="checkbox"/>	WEX5
I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,929.00 <input type="checkbox"/>	WEX5

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10076	Right Side Mount For Joystick	No Charge <input type="checkbox"/>	
I10077	Left Side Mount For Joystick	No Charge <input type="checkbox"/>	
I10074	Right Side Mount For Omni	No Charge <input type="checkbox"/>	
I10075	Left Side Mount For Omni	No Charge <input type="checkbox"/>	
I10078	Mount Seat Control Panel on Right	No Charge <input type="checkbox"/>	
I10079	Mount Seat Control Panel on Left	No Charge <input type="checkbox"/>	
I10080	Mount Seat Control Panel on Backrest	No Charge <input type="checkbox"/>	

## Alt Drive Systems Setup

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10576	(2) EXPANDABLE CONTROLLER (REQUIRED) Expandable Controller - Alternative Drive	No Charge <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10577	(3) ADJUSTABLE REMOVABLE HARDWARE (REQUIRED) Adjustable Removable Hardware - Alt Drive	110.00 <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
ID50881	(4) OMNI OR INPUT MODULE (REQUIRED) <i>Cannot select Input Module with Sip-n-Puff Input Devices (PN I10573 or I10574). Omni Control Module for R-net</i>	2,055.00 <input type="checkbox"/>	
I10581	<i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i> Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,450.00 <input type="checkbox"/>	

## Alt Drive Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
-------------	-------------	------------	----------

**Alt Drive Accessories**

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
IMICROSWITCH	Microlite Switch	148.00 <input type="checkbox"/>	
IBUDDYBUTTON	Buddy Button	111.00 <input type="checkbox"/>	
I5808	Mini Cup Switch	108.00 <input type="checkbox"/>	
IZTS-P	Proximity Switch	515.00 <input type="checkbox"/>	
IMINIJOY5SW	Mini Joystick 5 Switch OC	875.00 <input type="checkbox"/>	
I5018	Wafer Board w/5 Switches	497.00 <input type="checkbox"/>	

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code
- The client's measurements
- The client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil Ltd is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_