

# C350 PS



## Dealer Information

Contact: \_\_\_\_\_  
Dealer Code: \_\_\_\_\_  
Dealer Name: \_\_\_\_\_  
Address: \_\_\_\_\_  
City: \_\_\_\_\_  
Province/Postal Code: \_\_\_\_\_  
Phone#: \_\_\_\_\_  
Fax#: \_\_\_\_\_  
PO#: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*


First Name: \_\_\_\_\_  
Last Name: \_\_\_\_\_  
Diagnosis: \_\_\_\_\_  
Age: \_\_\_\_\_

## Client Measurements

**Weight:** \_\_\_\_\_

**Height:** \_\_\_\_\_

- A) Top of Shoulders: \_\_\_\_\_  
B) Chest Depth: \_\_\_\_\_  
C) Chest Width: \_\_\_\_\_  
D) Seat Depth: \_\_\_\_\_  
E) Top of Head: \_\_\_\_\_  
F) Elbow to Hand: \_\_\_\_\_  
G) Seat Pan to Elbow: \_\_\_\_\_  
H) Hip Width: \_\_\_\_\_  
I) Knee to Foot: \_\_\_\_\_



Display ADP Codes on order/quote: \_\_\_\_\_

Please send order/quote to fax#: **(800) 231-3256**

Canadian prices effective **October 28, 2008** and payable in Canadian dollars. All specifications and prices are subject to change without notice. **Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.**

**Permobil Ltd.**  
Toronto Canada  
C/O Permobil Inc  
6961 Eastgate Blvd.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
www.permobil.com

**Base***Select Options:*

Part Number	Description	MSRP - CAD	ADP Code
I10422	C350 Base 5mph - C350 PS1 <i>Price Includes: RWD C350 Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, VR2 Controller with 5mph Motor Programming, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Bodypoint Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,895.00 <input type="checkbox"/>	
I10554	C350 Base 6.5mph - C350 HS PS2 <i>Price Includes: RWD C350 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, 6.5mph Motor Programming, Tie-Down Hardware for Strap Systems, Anti-Tippers, Flat Free Tires, and Standard Bodypoint Positioning Belt. If selecting this Base, you MUST also select R-net Electronics. The 6.5mph Speed Package is NOT available with VR2 Electronics. Choose Flat Free or Pneumatic Tires Below.</i>	7,975.00 <input type="checkbox"/>	

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge <input type="checkbox"/>	
I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge <input type="checkbox"/>	

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I103532-99-0	C350 Shroud Vintage Wine	No Charge <input type="checkbox"/>	
I103533-99-0	C350 Shroud Crystal Blue	No Charge <input type="checkbox"/>	
I103529-99-0	C350 Shroud Onyx Black	No Charge <input type="checkbox"/>	
I103530-99-0	C350 Shroud Platinum Silver	No Charge <input type="checkbox"/>	
I103531-99-0	C350 Shroud Extreme Yellow	No Charge <input type="checkbox"/>	

*Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge <input type="checkbox"/>	
I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge <input type="checkbox"/>	

**Base Options***Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I103524-99-0	Power Adjustable Seat Height, 8" Travel <i>To assist in transfers and improve access for function.</i>	2,495.00 <input type="checkbox"/>	
I101994-99-0	Seat Tube Fixed	No Charge <input type="checkbox"/>	

*Select Options:*

Part Number	Description	MSRP - CAD	ADP Code
IM34 SLDG	Batteries - Group 34, 12V 60 Amp Hours, Installed	762.00 <input type="checkbox"/>	
I602307-99-0	P+ PP1A Programmer - VSI, VR2 or Pilot+ Controllers	595.00 <input type="checkbox"/>	
I103496-99-0	Light Kit for C350 <i>Cannot be ordered with VR2 Electronics.</i>	955.00 <input type="checkbox"/>	
ID51113	Programming Key for R-net - Permobil Version	435.00 <input type="checkbox"/>	

**Joystick Options***Must select one of the following:*

Part Number	Description	MSRP - CAD	ADP Code
I10045	Cone Shaped Joystick Knob - Standard.	No Charge <input type="checkbox"/>	
I100960-99-0	Chin Cup For Joystick Knob	94.00 <input type="checkbox"/>	

## Joystick Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I100961-99-0	"T" Handle For Joystick Knob	139.00 <input type="checkbox"/>	
I100962-99-0	"Large Ball" For Joystick Knob	94.00 <input type="checkbox"/>	
I100963-99-0	"Softball" For Joystick Knob	127.00 <input type="checkbox"/>	
I100964-99-0	"Mushroom" For Joystick Knob	127.00 <input type="checkbox"/>	
I100965-99-0	"Stick" For Joystick Knob	127.00 <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10424	Fixed Joystick Mount - Right - VR2	No Charge <input type="checkbox"/>	
I10425	Fixed Joystick Mount - Left VR2	No Charge <input type="checkbox"/>	
I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	325.00 <input type="checkbox"/>	WEPM
I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	325.00 <input type="checkbox"/>	WEPM
I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	375.00 <input type="checkbox"/>	WEPM
I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	375.00 <input type="checkbox"/>	WEPM
I101260-99-0	Midline Control Mount - Pilot+ & R-net Only	325.00 <input type="checkbox"/>	WE3D

Please indicate which armrest to mount midline below.

I10428	Mount Midline Joystick on Left Armrest	No Charge <input type="checkbox"/>	
I10427	Mount Midline Joystick on Right Armrest	No Charge <input type="checkbox"/>	
I10423	VR2 Basic Joystick - Standard	No Charge <input type="checkbox"/>	
I103544-99-0	R-net Remote Joystick w/ Mono Jacks <i>Expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	1,495.00 <input type="checkbox"/>	
ID51040	R-net Remote Joystick w/ Toggle Switches <i>Expandable/remote joystick with toggle switches. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	1,595.00 <input type="checkbox"/>	WEX1
ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle <i>Expandable/remote joystick with a single toggle switch for modes/profiles &amp; a speed pot. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	1,595.00 <input type="checkbox"/>	WEX1
I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i>	No Charge <input type="checkbox"/>	
ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 controllers.</i>	1,084.00 <input type="checkbox"/>	
I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	1,450.00 <input type="checkbox"/>	

## Seating System

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10281	PS Seating System Only - C300/C350 <i>Includes: Seat Frame, Fixed Back Frame, Fixed Legrests and Footplates, Fixed Arm Supports and Standard Armrest Pads. Weight Limit 300 lbs.</i>	No Charge <input type="checkbox"/>	
I10142	PS Seating System With 45 Degree CG Power Tilt <i>Price Includes: Seat Frame, All Seat Sizes, Angle Adjustable Back Canes, Height Adjustable Armrests, Standard Armrest Pads, Swing Away Legrests and Footplates. Weight Limit 300lbs.</i>	5,595.00 <input type="checkbox"/>	SEC1P0205

## Seating System Options

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10289	Single Seat Function Control Switch - LiberT <i>This device is a single toggle switch used to control only one seat function. If control of the seat function through the drive control is needed, then the drive electronics must be upgraded to Pilot+ or R-net Electronics, and a Single Seat Function Control Kit (Part Number I10129), must be selected.</i>	No Charge <input type="checkbox"/>	
I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	1,378.00 <input type="checkbox"/>	
I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	2,295.00 <input type="checkbox"/>	SEC1P1005
I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge <input type="checkbox"/>	
I102799-99-0	Toggle Switches on Seat Function Control Kit	245.00 <input type="checkbox"/>	
I103469-99-0	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge <input type="checkbox"/>	
I103470-99-0	Mount ICS Switchbox on Backrest	125.00 <input type="checkbox"/>	
I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,084.00 <input type="checkbox"/>	

## Backrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10144	Sling Upholstery Back Support - Standard <i>Width will be the same as seat width.</i>	No Charge <input type="checkbox"/>	
I10145	Curved Metal Back Support 22" High	410.00 <input type="checkbox"/>	
I10146	Curved Metal Back Support 22" High - No Charge <i>No charge if ordering tilt.</i>	No Charge <input type="checkbox"/>	
I10454	Omit Backrest	No Charge <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10147	Fixed Back Canes <i>Standard option.</i>	No Charge <input type="checkbox"/>	
I10148	Angle Adjustable Back Canes <i>Angle adjustable from 85 to 115 degrees.</i>	127.00 <input type="checkbox"/>	
I10149	Angle Adjustable Back Canes - No Charge <i>No charge if ordering tilt.</i>	No Charge <input type="checkbox"/>	

## Seat Sizing

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10150A	Seat Width 16"	No Charge <input type="checkbox"/>	
I10151A	Seat Width 18"	No Charge <input type="checkbox"/>	
I10219A	Seat Width 20"	No Charge <input type="checkbox"/>	
I10220A	Seat Width 22"	No Charge <input type="checkbox"/>	

### Seat Sizing

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10156	Seat Depth 16"	No Charge <input type="checkbox"/>	
I10157	Seat Depth 17"	No Charge <input type="checkbox"/>	
I10158	Seat Depth 18"	No Charge <input type="checkbox"/>	
I10159	Seat Depth 19"	No Charge <input type="checkbox"/>	
I10221	Seat Depth 20"	No Charge <input type="checkbox"/>	
I10222	Seat Depth 21"	No Charge <input type="checkbox"/>	
I10223	Seat Depth 22"	No Charge <input type="checkbox"/>	

### Armrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10164	Fixed Height Armrest Assembly PS <i>Flips back for transfers.</i>	No Charge <input type="checkbox"/>	
I10163	Height Adjustable Armrest Assembly PS <i>Adjustable in height, with a flip-back design.</i>	284.00 <input type="checkbox"/>	
I10230	Height Adjustable Armrest Assembly PS - No Charge <i>Adjustable in height, with a flip-back design. No charge if ordering power tilt.</i>	No Charge <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10290	PS Armrest Pads - 10.5"L x 2"W (Pair)	No Charge <input type="checkbox"/>	
I10291	PS Armrest Pads - 14"L x 2"W (Pair)	No Charge <input type="checkbox"/>	
I10292	Waterfall Armrest Pad for PS - 13.5"L x 3"W (Pair)	79.00 <input type="checkbox"/>	
I10168	Omit Armrest Pads	No Charge <input type="checkbox"/>	
I10166	Corpus Armrest Pad Adapter Brackets <i>Must be picked if ordering Corpus armrest pads.</i>	50.00 <input type="checkbox"/>	
I10208	Corpus Armrest Pad 4"x10" Left Leatherette	92.00 <input type="checkbox"/>	
I10209	Corpus Armrest Pad 4"x12" Left Leatherette	92.00 <input type="checkbox"/>	
I10210	Corpus Armrest Pad 4"x14" Left Leatherette	92.00 <input type="checkbox"/>	
I10211	Corpus Armrest Pad 4"x16" Left Leatherette	92.00 <input type="checkbox"/>	
I10212	Corpus Armrest Pad 4"x10" Right Leatherette	92.00 <input type="checkbox"/>	
I10213	Corpus Armrest Pad 4"x12" Right Leatherette	92.00 <input type="checkbox"/>	
I10214	Corpus Armrest Pad 4"x14" Right Leatherette	92.00 <input type="checkbox"/>	
I10215	Corpus Armrest Pad 4"x16" Right Leatherette	92.00 <input type="checkbox"/>	
I10169	Armrest Side Support Right Inner <i>Will limit the armrest from flipping back. For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	
I10170	Armrest Side Support Left Inner <i>Will limit the armrests from flipping back. For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	
I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	
I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	

### Legrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10172	70 Degree Swing Away Legrests (Pair)	199.00 <input type="checkbox"/>	
I10231	70 Degree Swing Away Legrests (Pair) <i>No charge if ordering Power Tilt.</i>	No Charge <input type="checkbox"/>	

## Legrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10316	Aeron Legrest Assembly for PS Seat <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes Corpus Style, Angle Adjustable, Two Piece Footplates.</i>	315.00 <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10173	Footplates (Pair), Standard Small	No Charge <input type="checkbox"/>	
I10174	Footplates (Pair), Standard Large <i>Not available with 16" seat width.</i>	No Charge <input type="checkbox"/>	
I10175	Small Footplates, Multi - Axis Adjustable, Aluminum	195.00 <input type="checkbox"/>	
I10176	Large Footplates, Multi - Axis Adjustable, Aluminum <i>Cannot be ordered with 16" or 18" seat widths.</i>	195.00 <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10177	Heel Loops	90.00 <input type="checkbox"/>	
I10197	Calf Supports Leatherette <i>Not available with lower leg lengths of 10" or less.</i>	No Charge <input type="checkbox"/>	
I10104	Calf Supports for Aeron Legrest Assembly - PS <i>Not available with lower leg lengths of 10" or less. Can be ordered only with Aeron Legrest Assembly for PS Seat (PN: I10316).</i>	No Charge <input type="checkbox"/>	

## Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10401	Bodypoint Non-Padded Hip Belt, PS - Medium	No Charge <input type="checkbox"/>	
I10406	Bodypoint Padded Hip Belt PS (C300) - Medium	119.00 <input type="checkbox"/>	
I10403	Bodypoint Non-Padded Hip Belt, PS - Large	119.00 <input type="checkbox"/>	
I10409	Bodypoint Padded Hip Belt PS - Large	129.00 <input type="checkbox"/>	
I10199	Headrest-Removable Leatherette <i>Must also select PN I10552.</i>	410.00 <input type="checkbox"/>	
I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	504.00 <input type="checkbox"/>	
I10180	Universal Headrest Adapter <i>Will be sent as part if ordering sling back upholstery.</i>	197.00 <input type="checkbox"/>	
I10552	Adj Removable - Headrest Hardware (1 Each)	110.00 <input type="checkbox"/>	SEHND0020
I102059-99-0	Medical Necessities Bag <i>Includes Back Pack Clips (Carabiner Style).</i>	265.00 <input type="checkbox"/>	
I102697001331	Medical Necessities Bag - Case Logic	265.00 <input type="checkbox"/>	
I10179	Back Pack Clips Carabiner Style <i>Will be sent as part if ordering sling back upholstery.</i>	75.00 <input type="checkbox"/>	
I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	356.00 <input type="checkbox"/>	SETND0010
I10556	Adj Removable - Tray Hardware (1 Each)	110.00 <input type="checkbox"/>	
I101964-99-0	Oxygen Holder for PS, Corpus or Standing Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc.</i>	224.00 <input type="checkbox"/>	WEPX
I101956-99-0	Crutch Holder for PS, Corpus or Standing Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder.</i>	224.00 <input type="checkbox"/>	

## Alternative Drive Systems

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	R-net Alternative Drive Input Devices		
	<i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Expandable Controller, (3) Adjustable/Removable Alt Drive Hardware, AND (4) Omni or Input Module.</i>		
	(1) INPUT DEVICES FOR ALTERNATIVE DRIVE CONTROL		
I10558	Gimble Joystick set up for Armrest Mount	1,388.00 <input type="checkbox"/>	WEX2
I10559	Gimble Joystick set up for Midline Mount	1,388.00 <input type="checkbox"/>	WEX2
I10560	Gimble Joystick set up for Swing Away Chin Boom	1,388.00 <input type="checkbox"/>	WEX2
I10561	<i>Custom charge may be applied if ordering lower backrest heights.</i> Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	1,555.00 <input type="checkbox"/>	WEX2
I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,555.00 <input type="checkbox"/>	WEX2
I10563	Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00 <input type="checkbox"/>	WEX2
I10564	Mini Joystick set up for Armrest Mount	2,035.00 <input type="checkbox"/>	WEX2
I10565	Mini Joystick set up for Midline Mount	2,035.00 <input type="checkbox"/>	WEX2
I10566	Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	2,035.00 <input type="checkbox"/>	WEX2
I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,135.00 <input type="checkbox"/>	WEX2
I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,135.00 <input type="checkbox"/>	WEX2
I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00 <input type="checkbox"/>	
I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,085.00 <input type="checkbox"/>	WEX3
I10570	Head Array - Adult 5 Switch HDR-A5	4,075.00 <input type="checkbox"/>	WEX4
I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	4,075.00 <input type="checkbox"/>	WEX4
I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	4,075.00 <input type="checkbox"/>	WEX4
I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,409.00 <input type="checkbox"/>	WEX5
I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,929.00 <input type="checkbox"/>	WEX5

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10076	Right Side Mount For Joystick	No Charge <input type="checkbox"/>	
I10077	Left Side Mount For Joystick	No Charge <input type="checkbox"/>	
I10074	Right Side Mount For Omni	No Charge <input type="checkbox"/>	
I10075	Left Side Mount For Omni	No Charge <input type="checkbox"/>	
I10078	Mount Seat Control Panel on Right	No Charge <input type="checkbox"/>	
I10079	Mount Seat Control Panel on Left	No Charge <input type="checkbox"/>	
I10080	Mount Seat Control Panel on Backrest	No Charge <input type="checkbox"/>	

## Alt Drive Systems Setup

**Alt Drive Systems Setup**

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	(2) EXPANDABLE CONTROLLER (REQUIRED)		
I10576	Expandable Controller - Alternative Drive	No Charge	<input type="checkbox"/>

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	(3) ADJUSTABLE REMOVABLE HARDWARE (REQUIRED)		
I10577	Adjustable Removable Hardware - Alt Drive	110.00	<input type="checkbox"/>

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	(4) OMNI OR INPUT MODULE (REQUIRED)		
	<i>Cannot select Input Module with Sip-n-Puff Input Devices (PN I10573 or I10574).</i>		
ID50881	Omni Control Module for R-net	2,055.00	<input type="checkbox"/>
	<i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>		
I10581	Input Module (IOM) for R-net Control System	1,450.00	<input type="checkbox"/>
	<i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>		

**Alt Drive Accessories**

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
IMICROSWITCH	Microlite Switch	148.00	<input type="checkbox"/>
IBUDDYBUTTON	Buddy Button	111.00	<input type="checkbox"/>
I5808	Mini Cup Switch	108.00	<input type="checkbox"/>
IZTS-P	Proximity Switch	515.00	<input type="checkbox"/>
IMINIJOY5SW	Mini Joystick 5 Switch OC	875.00	<input type="checkbox"/>
I5018	Wafer Board w/5 Switches	497.00	<input type="checkbox"/>

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code
- The client's measurements
- The client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil Ltd is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_