



Get Vertical

THE HELIUM is designed for active people for use indoors and outside. Constructed of high-strength aircraft aluminum and titanium, the Helium is the lightest weight stand-up wheelchair in the world. Brilliantly engineered, this ergonomically-friendly standing system incorporates gas struts that deliver the safe and easy lift you need.

A patented back angle compensation system and adjustable seat depth allow a perfect fit tailored to you. Sturdy and comfortable neoprene-foam knee blocks and an adjustable chest belt ensure the standing position is comfortable and secure.

The Helium increases independence and puts you eye-to-eye with family and friends. And it may provide a myriad of health benefits, including: better muscle tone and strengthened bones, improved circulation and digestion as well as reduced spasticity.



HELIUM



GET VERTICAL – Glide your way smoothly and easily to the standing position in a high-performance, light-weight wheelchair with clean lines and great looks to match.

LIFESTAND
BY PERMOBIL

HELIUM

STANDING INDEPENDENCE



Rear axle center of gravity adjustment for superior maneuverability



Durable, tension adjustable seat upholstery



Swing-away armrests (optional) allow obstacle-free transfers



Height and angle adjustable footplate for proper fit and comfort




Patented back angle compensation system allows proper positioning of the torso when standing



Simple disassembly for easy transporting



COLOR OPTIONS



White Chrome Titanium Anodized Blue Metal Black Metal Red

Seat Width	Seat Front to Floor Height	Rear Wheel Diameter	Front Wheel Diameter	Overall Width (0° camber)	Overall Width (3° camber)	Seat Depth (Knee front to backpost)	Seat Depth	Overall Length Seated	Overall Length Standing	Backrest Height	Overall Height (w/o push handles)	Max User Weight	Total Weight*	Weight of Heaviest Element
36cm/15"	51cm/20"	60cm/24"	12.5cm/5"	60cm/24"	66cm/26"	Adjustable 54cm/21" to 64cm/25"	Adjustable 43cm/17" to 54cm/21"	Adjustable 89cm/35" to 101cm/40"	Adjustable 93cm/37" to 105cm/41"	32cm/12.5" to 45cm/18"	78cm/31" to 91cm/36"	100 Kilos/ 220 lbs	From 16 Kilos/ 35 lbs	13.5 Kilos/ 30 lbs
38cm/15.5"	51cm/20"	60cm/24"	12.5cm/5"	62cm/24.5"	68cm/27"									
40cm/16"	51cm/20"	60cm/24"	12.5cm/5"	64cm/25"	70cm/27.5"									
42cm/16.5"	51cm/20"	60cm/24"	12.5cm/5"	66cm/26"	72cm/28.5"									
44cm/17.5"	51cm/20"	60cm/24"	12.5cm/5"	68cm/27"	74cm/29"									
46cm/18"	51cm/20"	60cm/24"	12.5cm/5"	70cm/27.5"	76cm/30"									

*dependent on sizes, types of wheels and options