



Pediatric Chair Descriptions



K450 MX

The all-wheel suspension system of the *new* K450 MX provides a smooth ride, whether cruising the halls at school or heading out to the playground. The power seat-to-floor function lets children get down on the ground, perfect for story time, transfers, and playing with others. The seat-to-floor height ranges from 3" to 26" making it easy to sit at a table, counter, and even reach things down low.

Plus, the K450 MX features R-net electronics

and has an array of power functions including: seat elevator, seat-to-floor, and optional power tilt. The modular seating system has center mount adjustable legrests and even comes with a free growth kit – something every kid needs.

K300 PS Jr.

The K300 comes with standard adult power base with suspension and the expandable seat pan allows for a custom fit no matter how fast a child grows. The low seat-to-floor height (16.5") allows children to interact with peers at their level. Features such as EZ adjust armrests, the power seat elevator and power tilt functions, R-Net electronics and the Uni-track system provide for maximum comfort and allow for an infinite number of adjustments and add-on accessories. The K300 offers easy touch swing-away footrests with adjustable high strength multi-axis footplates; available with 70° or 90° front hangers.



Koala Miniflex

Permobil's smallest pediatric wheelchair is designed with a small child in mind. It is an excellent choice to experience mobility for the first time. It is small, maneuvers easily and intuitively, and handles outdoors exceptionally well – all in all, it is a rugged little chair with the looks that are a perfect match for a child. The backward tilt is perfect when it is time for a little rest, or to realign position and get more comfortable. The five-degree forward tilt on the seat makes play easier and enhances alertness. The integrated, vertical height adjustment of the seat makes playing, bathing, grooming and eating easier. The Miniflex seat will

grow with the child for years to come. In addition to a variety of accessories, tubular interface is also available for use of aftermarket seating products.





Pediatric Chair Descriptions

C300 Corpus Jr.

The C300 Corpus Jr. offers unmatched growth capability and flexibility. The seat depth grows from 13" to 22", and the width can accommodate as narrow as 10" trunk and hip width yet grow all the way to 23". All seat functions are optional – it is available without functions, with 45° of tilt, recline, elevating legs or elevating seat. Any seat function can be added after the fact very quickly in the field. With tight turning, exceptional outdoors capabilities, flexibility of seat functions, ability to accept alternative drive controls and an economical base with low seat height, it is an ideal choice for a growing child.



C400 Corpus Jr.

The C400 Corpus Jr combines the flexibility of the Corpus Jr seating system with a compact, stable and sturdy base. Front wheel drive technology enables better outdoor performance while maintaining exceptional indoors maneuverability. Large batteries allow for prolonged driving before the risk of draining power, and four wheels enable better transfers. The Corpus is designed by one of the leading ergonomic experts in the world to match the natural contours of the body. The customization of the Corpus seating system is limitless. Depth, width, height, and angle of multiple seat components can be adjusted or customized. The seat depth grows from 13" to 22", and the width can accommodate as narrow as 10" trunk and hip width yet grow all the way to 23".



C400 VS Jr.

The *new* C400 VS Jr. redefines the way standing is incorporated into a child's daily routine and allows the user to drive while in a standing position. Standing can be achieved from either a sitting or a horizontal position, enabling the child to move about throughout the day and remain comfortable. Designed specifically for children, the C400 VS Jr. features a front-wheel drive base and a rugged suspension system. The small footprint of the C400 VS Jr. combines maneuverability with advanced obstacle climbing capabilities. Plus, the modular design of the Corpus seat ensures optimal fit with room for growth.

